

FULKU

FOREST UTILISATION LABORATORY, KYOTO UNIVERSITY

How to operate turnover_sensor in computer Clethra (R0.33)

Laboratory Analytical Procedure (LAP)

パソコンレスラでの turnover_sensor の操作方法 (R3)

実験室分析手順

Issue Date: January 2025

発行日：2025 年 1 月

Users: Eriko KITA, Yumi KIHARA, Shitephen WANG, Masako DANNOURA

FULKU is a laboratory based on Graduate school of Agriculture, Kyoto University.


Technical Report

Procedure Title: How to operate turnover_sensor in computer Clethra (R3)


Laboratory Analytical Procedure (LAP)

0. Move the “turnover_sensor.py” file in the root folder (Clethra)


(preset here in our lab → "C:\Users\Clethra")




1. open Anaconda Navigator




2. Chose the “arata” in Environments




3. Go “home” and launch CDM.exe Prompt



4. Type “python turnover_sensor.py” in cdm.exe and press enter



5. Select Folder (e.g. postprocess generated from ARATA)





6. Start to align images



6. Start to fitting images



7. Pick up the result in the pathway you set



Green is root growth; Red is root decomposed

Note:

1. If you have some error happened, it is good to ask ChatGPT first (e.g. the image sizes are different and the ARATA cannot output calculated values. ChatGPT will recommend you to remove images with different sizes)

2. If you want to know more details about the model, original source codes here <https://github.com/gn03138868/turnover-sensor>
3. If you still have some questions, you can ask Shitephen (gn03138868@gmail.com) via email
4. You can upgrade any packages and versions. It is just a small tool, easy to revise. (I developed it under the environment of ARATA.)

Python 3.7.16 (default, Jan 17 2023, 16:06:28) [MSC v.1916 64 bit (AMD64)]

Packages in python	Version
--------------------	---------

absl-py	1.4.0
anyio	3.5.0
argon2-cffi	21.3.0
argon2-cffi-bindings	21.2.0
astunparse	1.6.3
attrs	22.1.0
backcall	0.2.0
beautifulsoup4	4.11.1
bleach	4.1.0
brotlipy	0.7.0
cachetools	5.3.0
certifi	2022.12.7
cffi	1.15.1
charset-normalizer	3.1.0
colorama	0.4.6
cryptography	37.0.4
cuda-python	11.7.1
cvxopt	1.3.0
cycler	0.11.0
debugpy	1.5.1
decorator	5.1.1
defusedxml	0.7.1
entrypoints	0.4
ExifRead	3.0.0
fastjsonschema	2.16.2
flatbuffers	23.1.21
flit_core	3.6.0
fonttools	4.38.0
gast	0.4.0
google-auth	2.16.2
google-auth-oauthlib	0.4.6
google-pasta	0.2.0
grpcio	1.51.3
idna	3.4
importlib-metadata	6.0.0

importlib-resources	5.2.0
ipykernel	6.15.2
ipython	7.31.1
ipython-genutils	0.2.0
jedi	0.18.1
Jinja2	3.1.2
jsonschema	4.17.3
jupyter_client	7.4.9
jupyter_core	4.11.2
jupyter-server	1.23.4
jupyterlab-pygments	0.1.2
keras	2.11.0
kiwisolver	1.4.4
libclang	15.0.6.1
llvmlite	0.39.1
lxml	4.9.1
Markdown	3.4.1
MarkupSafe	2.1.1
matplotlib	3.5.3
matplotlib-inline	0.1.6
mistune	0.8.4
nbclassic	0.5.2
nbclient	0.5.13
nbconvert	6.5.4
nbformat	5.7.0
nest-asyncio	1.5.6
notebook	6.5.2
notebook_shim	0.2.2
numba	0.56.4
numpy	1.21.6
oauthlib	3.2.2
olefile	0.46
opt-einsum	3.3.0
packaging	22.0
pandas	1.3.5
pandocfilters	1.5.0
parso	0.8.3
patsy	0.5.3
pickleshare	0.7.5
Pillow	9.4.0
pip	22.3.1
pkgutil_resolve_name	1.3.10
prometheus-client	0.14.1
prompt-toolkit	3.0.36
protobuf	3.19.6

psutil	5.9.0
pyasn1	0.4.8
pyasn1-modules	0.2.8
pycparser	2.21
Pygments	2.11.2
pyOpenSSL	22.0.0
pyparsing	3.0.9
pyrsistent	0.18.0
PySocks	1.7.1
python-dateutil	2.8.2
pytz	2022.7.1
pywin32	305.1
pywinpty	2.0.10
pyzmq	23.2.0
requests	2.28.2
requests-oauthlib	1.3.1
rsa	4.9
scipy	1.7.3
Send2Trash	1.8.0
setuptools	65.6.3
six	1.16.0
sniffio	1.2.0
soupsieve	2.3.2.post1
statsmodels	0.13.5
tensorboard	2.11.2
tensorflow	2.11.0
tensorflow-estimator	2.11.0
tensorflow-intel	2.11.0
tensorflow-io-gcs-filesystem	0.31.0
termcolor	2.2.0
terminado	0.17.1
tinycss2	1.2.1
torch	1.13.1
torchaudio	0.12.0
torchvision	0.14.1
tornado	6.2
traitlets	5.7.1
typing_extensions	4.5.0
urllib3	1.26.15
wcwidth	0.2.5
webencodings	0.5.1
websocket-client	0.58.0
Werkzeug	2.2.3
wheel	0.38.4
win-inet-pton	1.1.0

wincertstore 0.2
zipp 3.11.0

5. Environment of User interface

Device name Clethra
Processor Intel(R) Core(TM) i5-9400F CPU @ 2.90GHz 2.90 GHz
Installed RAM 16.0 GB
Device ID B9F0C10B-4FA8-4D0B-9098-9A621D5DA7E8
Product ID 00325-81309-17231-AAOEM
System type 64-bit operating system, x64-based processor
Pen and touch No pen or touch input is available for this display

Edition Windows 11 Home
Version 22H2
Installed on 01/03/2023
OS build 22621.1555
Experience Windows Feature Experience Pack 1000.22640.1000.0